## **OtoAccess<sup>TM</sup>**

Interacoustics® Database



Supports all Interacoustics products




### **OtoAccess<sup>TM</sup>**


Interacoustics® Database

# User-friendly database

OtoAccess<sup>™</sup> is a user-friendly database for storage of audiologic data and client information in private and public ENT clinics. It is equally suitable for installation on PCs as a stand alone database application as it is for networks with a number of clients using the same SQL server.

Emphasis has been placed on high-speed data handling capacity. The number of patients and sessions is virtually unlimited and very high quantities of data will have little practical impact on search speed. OtoAccess™ has a journal feature for commenting on test sessions and a schedule planner for day-list management. Import and export facilities for XML and GDT formats are fully implemented. A high degree of data safety ensures that OtoAccess™ meets HIPAA requirements. Data transfer is possible from any version of IaBase to OtoAccess™.





#### **Additional features**

- High Speed Data Handling
- Unlimited patient storage
- Integrates multiple diagnostic procedures
- Journal entries for each session with standard word processing toolbar
- Import and Export Capabilities
- Print reports and EMR documents
- HIPAA Compliant
- Microsoft SQL Server 2005 express database for audiometric data storage
- Schedule planner with day-list
- 3 installation types: Stand-alone Network Client - Server/Network Client
- Data-transfer from any version of laBase to OtoAccess™





# Technical Specifications

| PC Minimum<br>requirements:  | <ul> <li>Processor of 800 MHz Pentium III or higher</li> <li>128 MB RAM or More (Windows® XP requires 128 MB RAM)</li> <li>Hard drive with 50MB free space for software installation and at least 44MB while running.</li> <li>Minimum display of 1024x768 pixels with 16bit color or more.</li> <li>One or more USB ports, version 1.1 or higher</li> <li>RS232 Serial port. (COM port for Audiometers -! Note Interacoustics can deliver USB to serial converters).</li> <li>Network LAN card (if connecting to database over the network)</li> <li>CD-ROM drive (convenience during installation and for CD writing back-up)</li> </ul> |
|------------------------------|--|
| Supported operating systems: | <ul> <li>Microsoft Windows® XP Professional (SP2 required)</li> <li>Microsoft Windows® XP Home Edition (SP2 required)</li> <li>Microsoft Windows® Server 2003 family</li> </ul>  |
| Software requirements:       | Microsoft Data Access Components 2.8   |
| Included parts:              | OtoAccess™ Database Software<br>Operation Manual on CD   |
| Optional parts:              | Diagnostic Suite Software  |

# All Interacoustics instruments supported:

- EP15 Clinical ABR
- EP25 Advanced ABR
- TEOAE25 Screening and diagnostic TEOAE
- DPOAE20 Screening and diagnostic DPOAE
- VEMP
- ABRIS Infant screening
- Affinity Hearing Aid Analyzer
- Audiometers
- Middle Ear Analyzers



300/302 - 2 - 10/200

#### Sales and service in your area:

Read more here: www.interacoustics.com/com/OtoAccess



Phone: +45 6371 3555 · Fax: +45 6371 3522

E-mail: info@interacoustics.com DK-5610 Assens, Denmark Web: www.interacoustics.com